<b>ico   cet</b>	<b>(Frépaño)</b> BAAT Medical Products B.V., F. Haezeleijstraat 800, 7555 RJ Hengelo, NL Telefónu: +31-88-56 56 600 www.baaomedical.com <b>Distribuidor:</b> icotec ag, Industriestrasse 12, 9450 Abtstetten, Switserlan Phone: +41 71 757 00 00, info@icotec.ch, www.icotec-medical.com
180198190710_FU_rev2(2020)	

<b>ENGLISH</b>
<b>IFU</b> <b>Trépaño (BAAT Medical)</b> for removal of icotec Carbon/PEEK pedicle screws
<b>INTENDED USE</b> The Trepanes are intended for the removal of icotec Carbon/PEEK pedicle screws when the sticks or broken and cannot be extracted with standard instrumentation.
<b>CAUTION:</b> Federal (USA) law restricts this device to use by the order of a physician.
<b>USE INFORMATION</b> The Trepanes must only be used by surgeons who are fully familiar with the surgical technique required. The operating surgeon must take care not to use the instruments to exert inappropriate stress on the patient or the implants and must scrupulously comply with any operating procedure described in the surgical technique provided. For example, the forces exerted when repositioning an instrument in-site must not be excessive as this is likely to cause injury to the patient.
To reduce the risks of breakage, it is not allowed to impact the instrument. Extreme care must be taken when the devices are used near vital organs, nerves or vessels.
The surgical steps and warnings can be found in the BAAT Medical Surgical Technique (EN_180198190710_ST_rev1(2020)).
<b>LIMITATIONS ON REPROCESSING</b> The Trepanes are intended for use on a single patient during a single procedure. Once used, they must be discarded. Reprocessing is not allowed if the Trepane was used.

<b>CLEANING AND DISINFECTION</b> <b>The Trepanes are made of stainless steel (1, 4020)</b>
<b>Manual pre-cleaning</b>
<ul style="list-style-type: none"><li>Wipe under cool tap water (approx. 18 °C) for 10 s</li> <li>Sonicate in cleaning solution 0.5% Neodisher MedClean Forte (Dr. Weigert) at 40 °C for 5 min</li> <li>Treat the outer surface under cool tap water with a nylon brush until visibly clean</li> <li>Treat the inner surface/lumina under cool tap water with a white tooth brush until visibly clean</li> <li>Rinse with cold desalinated water for 10 s</li></ul>
<b>Automated cleaning and disinfection</b>
Use a washer-disinfector according ISO 15883-1/2. Load products as follows. Use the following information to guarantee an optimal rinsing. Use the camulations and holes can drain. The trays must not be overloaded to guarantee an optimal rinsing. Use the instructions in the manual of the washer-disinfector. <ul style="list-style-type: none"><li>2 min pre-cleaning with cool tap water</li> <li>5 min cleaning with 55 °C tap water and 0.5% Neodisher MedClean Forte</li> <li>Draining</li> <li>3 min rinsing with cold desalinated water</li> <li>Draining</li> <li>2 min rinsing with cold desalinated water</li> <li>Draining</li> <li>5 min thermal disinfection with 90 °C desalinated water (pH-value &gt; 3000).</li></ul> <p>Other parameters are acceptable when an AD-value of 3000 is achieved.</p> <ul style="list-style-type: none"><li>Draining</li> <li>Drying, do not exceed 120 °C</li></ul>
<b>MAINTENANCE AND INSPECTION</b>
The product must be examined for visible damage such as cracks, deformations, wear and corrosion. Cutting edges should be free of nicks and present a continuous edge. Discard blunt or damaged instruments.

<b>STERILIZATION</b>
Products may be loaded into the dedicated tray of a general purpose sterilization tray • Disposable area in sterilization paper according ISO 11607-1, EN 868-2
<b>Steam sterilization</b>
When sterilizing multiple products in one autoclave cycle ensure that the sterilizer's maximum load is not exceeded. <ul style="list-style-type: none"><li>Method: Pre-vacuum dynamic air-removal according EN 13060 / EN 285 / ISO 17865</li> <li>Exposure time: 4 minutes or</li> <li>Method: Pre-vacuum dynamic air-removal according EN 13060 / EN 285 / ISO 17865</li> <li>Temperature: 132 °C</li> <li>Exposure time: 3 minutes</li></ul>
<b>Drying</b>
The integrity of packaging and containers should be visually checked after removal from the sterilizer. Damaged packaging and containers should be treated as non-forming product. Drying should be carried out in an environment in which particles and microbial contamination are controlled. <ul style="list-style-type: none"><li>Drying time: 20 minutes</li></ul>
<b>STORAGE</b>
Products must be stored in a controlled environment.
<b>DISPOSAL</b>
The disposal of this medical product requires no special measures. Be sure to observe all national/local regulations and guidelines when disposing of the packaging material and potentially infectious items.

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